



# AquaFlow™


Bubble humidifier range



Oxygen and Aerosol Therapy ▪ Humidification

## AquaFlow™ oxygen bubble humidifiers

A comprehensive range of devices for delivering unheated humidification to your patients. Supplied with a bottle, the AquaFlow range of bubble humidifiers are suitable for patients receiving variable oxygen therapy.



| Code    | Description   | Box Qty. |
|---------|---|----------|
| 1507000 | AquaFlow, oxygen bubble humidifier with bottle, 500 ml                    | 20       |
| 1505000 | AquaFlow, oxygen bubble humidifier with bottle plus M12 adaptor, 500 ml   | 20       |
| 1521000 | AquaFlow, homecare oxygen bubble humidifier with bottle, 120 ml           | 20       |
| 1520000 | AquaFlow, homecare bubble humidifier with bottle plus M12 adaptor, 120 ml | 20       |

Make an enquiry

## Flowmeter adaptors

Oxygen therapy flowmeters are available with a variety of connection types. The Intersurgical range of adaptors has been designed to permit the attachment of humidification devices with 9/16 UNF Diameter Index Safety System (D.I.S.S) threaded couplings to a variety of flowmeter types.

Please download the explanatory poster entitled Oxygen Therapy Flowmeters to choose the most suitable adaptor for your application:



[www.intersurgical.com/products/oxygen-and-aerosol-therapy/flowmeter-adaptors](http://www.intersurgical.com/products/oxygen-and-aerosol-therapy/flowmeter-adaptors)

| Code    | Description            | Box Qty. |
|---------|------------------------|----------|
| 1506001 | BOC® flowmeter adaptor | 50       |
| 1509001 | Flowmeter adaptor M12  | 50       |

Make an enquiry

BOC is a registered trademark of The BOC Group Ltd.

IS10.10 • Issue 12 09.24



Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK  
T: +44 (0)118 965 6300 info@intersurgical.com www.intersurgical.com



The manufacturer Intersurgical Ltd is certified to ISO 14001:2015, ISO 9001:2015, ISO 13485:2016 and MDSAP

**Please think before you print**  
Save energy and paper.  
If you must print this information sheet please print it double sided.